

Statement regarding the NucleoCounter system given by:

Corey W. Pine
Research Associate
Amicus Therapeutics, Inc.
6 Cedarbrook Dr
Cranbury, NJ 08512
PH : 609-662-2024

“The NucleoCounter is a sleek and convenient alternative to the subjective method of the hemacytometer. The ease of use and lack of cleaning or maintenance makes cell biology considerably easier and more efficient. I have used the competition's particle counters (Coulter) and have found the cleaning and maintenance almost as time consuming as counting by hand. Another advantage of the PI based counting system over the various particle counters is the inability of the trypan-blue based systems to differentiate cells from air bubbles and accurately count aggregated or agglutinated cell clusters. I am more than pleased with my company's purchase of the NucleoCounter and have recommended it to my colleagues in the academic community.”

Corey W. Pine